

17 June 2022

New Trace Element EDTA tubes

From late June, SA Pathology will stock BD Trace Element EDTA tubes (navy top) for Zinc, Copper, Selenium, Cobalt and Chromium analyses.



Oracle order number: 19050009

The new tube provides improved detection of trace element levels and reduced risk of false high Zinc results.

The updated Order of Draw-Quick Guide (PUB-0085) is attached for your reference and is also available on the SA Pathology <u>website</u>.

Note – this change does not apply to paediatric collections.

For more information, please contact Chemical Pathology on 08 8222 3107.

Regards,

Dr Sophia OttoClinical Service Director
SA Pathology

www.sapathology.sa.gov.au

